




Ventilators on the Fly


4 Basic Questions:

What is the mode?

What is the level of oxygen support?

Is ventilation appropriate?

How are the respiratory mechanics?



Authors: David Scudder, MD; Kaitlyn McLeod, MD; Arun Kannappan, MD

Section Editor: Trevor Steinbach, MD

Executive Editor: Yilin Zhang, MD

Orientation

Mode


O₂ support

Ventilation

Mechanics

Cases

VENTILATOR ORIENTATION



Actual values measured
by the ventilator

Set by the
clinician

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

General Approach

- 1 Classic square wave
- 2 Mandatory RR?
- 3 What you target
- 4 What you monitor
- 5 When to use


WHAT IS THE MODE?

Pressure Control

Volume Control

PRVC

Pressure Support



Click for summary table of ventilator modes

- Orientation
- Mode
- O₂ support
- Ventilation
- Mechanics
- Cases

WHAT IS THE MODE?


- Pressure Control
- Volume Control
- PRVC
- Pressure Support

Key outputs are tidal volume and minute ventilation

Pressure Control Mode

- Classic square wave
- Mandatory RR?
- What you target
- What you monitor
- When to use

Patient comfort, precise airway pressure adjustment



Summary of ventilator modes

WHAT IS THE MODE?

Pressure Control

Volume Control

PRVC

Pressure Support

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

Key output is
peak pressure

Volume Control Mode

Classic square wave

Mandatory RR?

What you target

What you monitor

When to use



ARDS (lung-protective tidal volumes)

P_{Peak}
22 cm H₂O

V_T
360 mL

RR
25 min⁻¹

Min Vent
9.0 L/min



24

RR

"Control
mode" =
mandatory
RR

360

V_T

Set target
tidal volume
(V_T)

5

PEEP

40

FiO₂

Summary of ventilator modes

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

Pressure Regulated Volume Control (PRVC)

Combined mode. Clinician sets a target V_T . Ventilator uses an algorithm to determine the pressure required to achieve target.

What you monitor

When to use

Lung protective tidal volumes with
comfort of pressure support

WHAT IS THE MODE?

Pressure Control



Volume Control

PRVC

Pressure Support

Key outputs
are peak
pressure and
achieved tidal
volume

P_{Peak}
22 cm H₂O
 V_T
358 mL
RR
25 min⁻¹
Min Vent
9.0 L/min



24

RR

“Control
mode” =
mandatory
RR

360

V_T

Set target
tidal volume
(V_T)

5

PEEP

40

FiO₂

Summary of ventilator modes

- Orientation
- Mode
- O₂ support
- Ventilation
- Mechanics
- Cases

WHAT IS THE MODE?


- Pressure Control
- Volume Control
- PRVC
- Pressure Support

Key outputs are tidal volume, RR, and minute ventilation

Pressure Support Mode

- Classic square wave
- Mandatory RR?
- What you target
- What you monitor
- When to use

SBTs, patient providing adequate RR



Summary of ventilator modes

VENTILATOR MODE SUMMARY

Pressure Control

Volume Control

PRVC

Pressure Support

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

Mode	Classic square wave	Mandatory RR?	What you target	What you monitor	When to use
Pressure Control	Pressure	Yes	Inspiratory Pressure (P_{Control})	Tidal Volume (V_T) and Minute Ventilation (V_E)	Patient comfort, precise airway pressure adjustment
Volume Control	Flow	Yes	Tidal Volume (V_T)	Peak Pressure and Plateau Pressure	ARDS (lung-protective tidal volumes)
Pressure Regulated Volume Control	Pressure	Yes	Tidal Volume (V_T)	Peak Pressure, Plateau Pressure, Achieved Tidal Volume (V_T)	Lung protective tidal volumes with comfort of pressure support
Pressure Support	Pressure	No	Pressure support (P_{Support})	Tidal Volume (V_T), Respiratory Rate (RR), and Minute Ventilation (V_E)	SBTs, patient providing adequate RR

WHAT IS THE LEVEL OF O₂ SUPPORT?



Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

FiO ₂	PEEP
30%	5
50%	8
70%	10
100%	18-24

P_{Peak}
22 cm H₂O
V_T
355 mL
RR
25 min⁻¹
Min Vent
8.9 L/min

Pressure

Flow

Time

Time

24

RR

360

V_T

5

PEEP

40

FiO₂

Higher PEEP
and FiO₂ =
higher O₂
support

The ARDSNet PEEP ladder gives a sense of where the patient is with oxygen support.

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases


IS VENTILATION APPROPRIATE?

1. Check the minute ventilation

2. Is it appropriate?

3. Vent or Patient?

P_{Peak}
22 cm H₂O
 V_T
355 mL
RR
20 min⁻¹
Min Vent
7.1 L/min



18

RR

360

V_T

5

PEEP

35

FiO₂

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

IS VENTILATION APPROPRIATE?

1. Check the minute ventilation


2. Is it appropriate?

3. Vent or Patient?

$$\text{Min Vent (V}_E\text{)} = V_T \times \text{RR}$$

$$7.1 \text{ L/min} = 0.355 \text{ L} \times 20 \text{ min}^{-1}$$

Normal is 4-6 L/min
Often much higher when
disease is present



18

RR

360

V_T

5

PEEP

35

FiO₂

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

pH/pCO₂/pO₂/HCO₃


General rule is to maintain
pH 7.2-7.5

IS VENTILATION APPROPRIATE?

1. Check the minute ventilation

2. Is it appropriate?

3. Vent or Patient?



Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

IS VENTILATION APPROPRIATE?

1. Check the minute ventilation

2. Is it appropriate?

3. Vent or Patient?

Is patient triggering breaths?

no

Adjust ventilator settings


$$V_E = V_T \times RR$$

yes

Look for patient causes

CNS suppression
CNS injury
Pain
Salicylates

P_{Peak}
22 cm H₂O
V_T
355 mL
RR
20 min⁻¹
Min Vent
7.1 L/min



18

RR

360

V_T

5

PEEP

35

FiO₂

- Orientation
- Mode
- O₂ support
- Ventilation
- Mechanics**
- Cases

HOW ARE THE RESPIRATORY MECHANICS?

Pressure-Targeted Modes

Key Parameters: V_T and V_E

Low V_T and V_E point to trouble with respiratory mechanics

P_{Peak}
28 cm H₂O

V_T
125 mL


RR
25 min⁻¹

Min Vent
3.2 L/min

Volume-Targeted Modes

Key Parameter: ΔP

Most lungs with good mechanics require a ΔP of <10 cm H₂O to achieve a 6 mL/kg tidal volume



What if $\Delta P > 10$?

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

HOW ARE THE RESPIRATORY MECHANICS?


In Volume-targeted modes, if $\Delta P > 10$, perform an inspiratory hold

P_{Peak}
29 cm H₂O

V_T
355 mL

RR
24 min⁻¹

Min Vent
8.5 L/min



24

RR

360

V_T

8


PEEP

35

FiO₂

HOW ARE THE RESPIRATORY MECHANICS?

In Volume-targeted modes, if $\Delta P > 10$, perform an inspiratory hold



Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

Inspiratory hold

$P_{Plateau}$

Resistance Issue


Compliance Issue

$Peak - Plat > 3$

Mucous plug
Kinked tube
Bronchospasm


$Plat - PEEP > 10$

Pneumothorax
Alveolar filling
Fibrosis
Mainstem intubation
Distended abdomen



Case 1

- What is the mode?
Pressure control
- What is target parameter?
P_{Control}
- What are the key outputs?
V_T and minute ventilation
- What is the actual RR?
21 breaths per minute
- How much O₂ support is the patient on?
Low O₂ support – patient is on minimal FiO₂ and PEEP



Case 2

Orientation

Mode

O₂ support

Ventilation

Mechanics

Cases

1

2

3

What is the mode?

Volume control

What is target parameter?

V_T

What are the key outputs?

Peak pressure ($\Delta P < 10$, so
inspiratory hold not needed)

What is the actual RR?

25 breaths per minute

How much O₂ support is the
patient on?

Very high O₂ support.
Maximum FiO₂ and PEEP.

4

P_{Peak}
27 cm H₂O

V_T
355 mL

RR
25 min⁻¹

Min Vent
8.9 L/min

Pressure

$\Delta P = \text{Peak} - \text{PEEP}$

Time

Flow

1

Square wave in flow =
volume mode

Time

2

24

RR

3

360


V_T

18

PEEP


100

FiO₂



- What is the mode?
Pressure support
- What is target parameter?
P_{Support}
- What are the key outputs?
V_T, RR, and minute ventilation
- What is the actual RR?
25 breaths per minute
- How much O₂ support is the patient on?
Low O₂ support. Low FiO₂ and PEEP

Case 3



Modern Framing — Lung-Protective Ventilation (unchanged core, 2015–24 evidence)

The original talk teaches mode identification; this slide adds current ARDS ventilation targets. The standing teaching is unchanged.

Low tidal volume (still the foundation)

- 6 mL/kg predicted body weight (range 4–8) in ARDS; the original ARDSNet mortality benefit (NEJM 2000) remains standard of care.

Pressure limits

- Keep plateau pressure <30 cm H₂O (many target ≤28).

Driving pressure (ΔP = Plateau – PEEP) — the newer prognostic lever

- Lower driving pressure is the ventilator variable most strongly associated with survival in ARDS (Amato, NEJM 2015; mediation analysis of 3,562 patients). Target ΔP ≤14–15 cm H₂O where achievable.
- Note: on this slide ΔP = Peak – PEEP (a quick bedside screen). True driving pressure uses the PLATEAU pressure (inspiratory hold), not peak.

Sources: ARDSNet, NEJM 2000;342:1301. Amato MBP et al, NEJM 2015;372:747. Lung-protective targets unchanged through 2024.

Take-Home Points

- Read any ventilator with 4 questions: (1) What mode? (2) What O₂ support? (3) Is ventilation appropriate? (4) How are the mechanics?
- Mode clues: square wave in PRESSURE = pressure mode (PC, PRVC, PS); square wave in FLOW = classic volume control. A set/mandatory RR = control mode; no mandatory RR = support mode (PS).
- Target vs monitor: set pressure → monitor VT & minute ventilation; set volume → monitor peak & plateau pressure. PEEP and FiO₂ set the O₂ support in every mode (ARDSNet PEEP/FiO₂ ladder).
- Ventilation: $VE = VT \times RR$ (normal 4–6 L/min, often higher in disease); aim to keep pH ~7.2–7.5. Check whether breaths are vent- or patient-triggered.
- Mechanics: in volume-targeted modes, if $\Delta P > 10$ do an inspiratory hold. Peak – Plateau > 3–5 = RESISTANCE problem (mucus plug, kinked tube, bronchospasm); elevated Plateau – PEEP (>10) = COMPLIANCE problem (pneumothorax, edema/alveolar filling, fibrosis, mainstem intubation, distended abdomen).

Attribution & Changelog

Original

- "Ventilators on the Fly" — Scudder D, McLeod K, Kannappan A; sec. ed. Steinbach T; exec. ed. Zhang Y. TeachIM.org, Dec 2021.
https://teachim.org/teaching_material/ventilators-on-the-fly/ — CC BY-NC 4.0. Educational use only; not medical advice.

Correction made to this copy (June 2026)

- Slide 16: "Compliance Issue" box label "Plat – Peak > 10" → "Plat – PEEP > 10" — fixes an internal inconsistency with the same slide's waveform annotation; compliance is flagged by elevated plateau relative to PEEP.

Appended (originals untouched)

- Lung-protective / driving-pressure framing slide; take-home points; this attribution slide. All interactive animations and click-to-reveal builds preserved.

Companion files

- Learner Handout (PDF) and Teachers Guide (DOCX) copied unchanged into the lesson folder & linked — reviewed, no factual errors.